


CONTENTS

<i>Needle Holders</i>	pages	2-12
<i>Forceps</i>	pages	13-37
<i>dressing</i>	pages	14-16
<i>tissue</i>	pages	16-20
<i>atraumatic</i>	pages	20-25
<i>tying</i>	pages	26-31
<i>jewelers</i>	pages	31-33
<i>different</i>	pages	34-37
<i>Scissors</i>	pages	38-48
<i>Various instruments</i>	pages	49-55
<i>hooks</i>	pages	50
<i>retractors</i>	pages	51,55
<i>probes</i>	pages	52,54
<i>loops</i>	pages	53
<i>dissectors</i>	pages	53,54
<i>vascular clips</i>	pages	55
<i>bulldog clamps</i>	pages	55
<i>Numerical index</i>	page	56

TITANIUM NEEDLE HOLDERS



If you want to order the Needle Holder with:
a) regular lock - add "L" to the instrument code | VH 600 L
b) universal lock - add "U" to the instrument code | VH 600 U





TC-Tungsten Carbide

Tungsten carbide particles metallurgically fused to the jaw create a fine textured surface.


This extremely durable surface allows for the secure, atraumatic handling of delicate suture and tissue.

NEEDLE HOLDERS


VH 600




VH 600.01




VH 600.02



VH 600.03





VH 600.04




Tungsten carbide jaws

Overall length
135 mm 5.3 inches
Round handle d=9 mm
Universal lock

VH 601



VH 601.01




Tungsten carbide jaws

Overall length
135mm 5.3 inches
Round handle d=6 mm
Universal lock



TITANIUM

TITANIUM

VH 602




VH 602.01




T u n g s t e n c a r b i d e j a w s

Needle Holder-Scissors
Overall length
140 mm 5.5 inches
Round handle d=8 mm

VH 603



VH 603.01



T u n g s t e n c a r b i d e j a w s



Overall length
145 mm 5.7 inches
Round handle d=8 mm
Universal lock

TITANIUM


TITANIUM

NEEDLE HOLDERS


VH 604




VH 604.01



VH 604.02





VH 604.03




Tungsten carbide jaws

Overall length
150 mm 5.9 inches
Round handle d=10 mm
Universal lock

VH 605



VH 605.01



Tungsten carbide jaws



Overall length
160 mm 6.3 inches
Round handle d=9 mm
Universal lock

TITANIUM


TITANIUM

NEEDLE HOLDERS


VH 606




VH 606.01




VH 606.02




VH 606.03



VH 606.04





VH 606.05




Tungsten carbide jaws

Overall length
165 mm 6.5 inches
Flat handle 8 mm wide
Universal lock


VH 607




VH 607.01




VH 607.02




VH 607.03



VH 607.04



VH 607.05



Tungsten carbide jaws


Overall length
165 mm 6.5 inches
Round handle d=9 mm
Universal lock

TITANIUM



TITANIUM

NEEDLE HOLDERS

VH 608




VH 608.01




Tungsten carbide jaws

Overall length
180 mm 7.1 inches
Round handle d=9 mm
Universal lock


VH 609




VH 609.01




VH 609.02




VH 609.03




VH 609.04





VH 609.05



VH 609.06



VH 609.07





Tungsten carbide jaws

Overall length
185 mm 7.3 inches
Round handle d=8 mm
Universal lock


TITANIUM

TITANIUM


VH 610




VH 610.01




VH 610.02




VH 610.03



VH 610.04





VH 610.05




Tungsten carbide jaws

Overall length
185 mm 7.3 inches
Flat handle 8 mm wide
Universal lock

VH 611



VH 611.01



Tungsten carbide jaws


Overall length
185 mm 7.3 inches
Round handle d=9 mm
Universal lock

TITANIUM


TITANIUM

NEEDLE HOLDERS


VH 612



VH 612.01




T u n g s t e n c a r b i d e j a w s




Overall length
200 mm 7.9 inches
Round handle d=9 mm
Universal lock


VH 613




VH 613.01




VH 613.02



VH 613.03



T u n g s t e n c a r b i d e j a w s




Overall length
200 mm 7.9 inches
Round handle d=9 mm
Universal lock


TITANIUM

TITANIUM


VH 614




VH 614.01



VH 614.02




Tungsten carbide jaws




Overall length
200 mm 7.9 inches
Flat handle 12 mm wide


VH 615




VH 615.01




VH 615.02



VH 615.03



Tungsten carbide jaws





Overall length
210 mm 8.3 inches
Round handle d=9 mm
Universal lock

TITANIUM


TITANIUM

NEEDLE HOLDERS


VH 616




VH 616.01




VH 616.02




VH 616.03



VH 616.04





VH 616.05




Tungsten carbide jaws

Overall length
210 mm 8.3 inches
Flat handle 9 mm wide
Universal lock


VH 617




VH 617.01




VH 617.02




VH 617.03



VH 617.04



VH 617.05




Tungsten carbide jaws

Overall length
210 mm 8.3 inches
Round handle d=8 mm
Universal lock


TITANIUM

TITANIUM


VH 618




VH 618.01




VH 618.02




VH 618.03



VH 618.04




VH 618.05




Tungsten carbide jaws

Overall length
230 mm 9.1 inches
Round handle d=9 mm
Universal lock


VH 619




VH 619.01




VH 619.02




VH 619.03



VH 619.04



VH 619.05



Tungsten carbide jaws

Overall length
240 mm 9.4 inches
Flat handle 9 mm wide
Universal lock

TITANIUM

TITANIUM

TITANIUM FORCEPS



dressing


tissue

atraumatic

tying

jewelers

different




TC-Tungsten Carbide

Tungsten carbide particles metallurgically fused to the jaw create a fine textured surface.



This extremely durable surface allows for the secure, atraumatic handling of delicate suture and tissue.

Dressing

VF 700




VF 700.01




Adson
Overall length
120 mm 4.7 inches
Flat handle 14 mm wide

Dressing


VF 701




VF 701.01



VF 701.02




VF 701.03



Adson
Overall length
130 mm 5.1 inches
Flat handle 14 mm wide

Dressing

VF 702




T u n g s t e n c a r b i d e j a w s


Sharp Pointed Tips
Overall length
145 mm 5.7 inches
Flat handle 14 mm wide

Dressing

VF 703




VF 703.01




Serrated platform
Overall length
152 mm 5.9 inches
Flat handle 7.4 mm wide

Dressing


VF 704




VF 704.01



VF 704.02




VF 704.03



Overall length
200 mm 7.8 inches
Round handle d=9 mm

Tissue

VF 705



Tungsten carbide jaws

Colibri Style
Overall length
120 mm 4.7 inches
Round handle d=8 mm


TITANIUM

TITANIUM


FORCEPS

Tissue


VF 706




VF 706.01



VF 706.02




Tungsten carbide jaws




Overall length
120 mm 4.7 inches
Flat handle 14 mm wide

Tissue

VF 707



Tungsten carbide jaws




Overall length
130 mm 5.1 inches
Round handle d=8 mm

TITANIUM


TITANIUM

Tissue

VF 708



VF 708.01




Tungsten carbide jaws

Overall length
130 mm 5.1 inches
Flat handle 14 mm wide

Tissue

VF 709



Tungsten carbide jaws

Overall length
145 mm 5.7 inches
Flat handle 14 mm wide



TITANIUM

TITANIUM

FORCEPS

Tissue

VF 710





T u n g s t e n c a r b i d e j a w s

Overall length
180 mm 7.0 inches
Round handle d=8 mm


Tissue

VF 711



T u n g s t e n c a r b i d e j a w s

VF 711.01





Overall length
200 mm 7.8 inches
Round handle d=9 mm

TITANIUM

TITANIUM

Tissue

VF 712




T u n g s t e n c a r b i d e j a w s

Bayonet
Overall length
200 mm 7.8 inches
Flat handle 12 mm wide

Atraumatic

VF 713



Overall length
130 mm 5.7 inches
Flat handle 14 mm wide



TITANIUM

TITANIUM

FORCEPS

Atraumatic


VF 714



De Bakery
Overall length
150 mm 5.9 inches
Round handle d=8 mm


Atraumatic

VF 715




Tungsten carbide jaws


VF 715.01




VF 715.02



VF 715.03




VF 715.04



De Bakery
Overall length
150 mm 5.9 inches
Flat handle 11 mm wide


Atraumatic

VF 716




2.5x2.0 mm

VF 716.01



1.5x1.8 mm

VF 716.02




2.5x2.0 mm




De Bakey
Overall length
185 mm 7.2 inches
Flat handle 11 mm wide

Atraumatic

VF 717



1.0x1.0 mm



De Bakey
Overall length
185 mm 7.2 inches
Round handle d=9 mm


TITANIUM

TITANIUM


FORCEPS

Atraumatic

VF 718



VF 718.01





De Bakey
Overall length
200 mm 7.8 inches
Flat handle 11 mm wide



Atraumatic

VF 719



De Bakey
Overall length
200 mm 7.8 inches
Round handle d=9 mm





TITANIUM

TITANIUM

Atraumatic



VF 720



De Bakey
Overall length
210 mm 8.2 inches
Round handle d=9 mm

Atraumatic

VF 721



De Bakey (Bayonet)
Overall length
200 mm 7.8 inches
Flat handle 12 mm wide



TITANIUM

TITANIUM

FORCEPS

Atraumatic


VF 722




De Bakey (Bayonet)
Overall length
220 mm 8.6 inches
Flat handle 12 mm wide

Atraumatic


VF 723




VF 723.01



VF 723.02



VF 723.03





De Bakey
Overall length
240 mm 9.4 inches
Flat handle 13 mm wide

TITANIUM

TITANIUM

Tying

VF 724





Tungsten carbide jaws

Sharp Pointed Tips
Overall length
120 mm 4.7 inches
Round handle d=8 mm


Tying

VF 725




Tungsten carbide jaws


VF 725.01



VF 725.02



VF 725.03




Overall length
135 mm 5.3 inches
Round handle d=8 mm


FORCEPS

Tying


VF 726



VF 726.01




Tungsten carbide jaws




Overall length
152 mm 5.9 inches
Round handle d=8 mm

Tying


VF 727




VF 727.01




VF 727.02



VF 727.03



Tungsten carbide jaws




Overall length
155 mm 6.1 inches
Round handle d=8 mm

TITANIUM


TITANIUM

Tying


VF 728




VF 728.01




VF 728.02



VF 728.03



VF 728.04




Tungsten carbide jaws


Overall length
160 mm 6.2 inches
Round handle d=8 mm

Tying


VF 729




VF 729.01



VF 729.02



VF 729.03




Tungsten carbide jaws

Overall length
175 mm 6.8 inches
Round handle d=8 mm


FORCEPS

Tying


VF 730



VF 730.01




Tungsten carbide jaws




Overall length
178 mm 7.0 inches
Round handle d=8 mm

Tying


VF 731



VF 731.01



Tungsten carbide jaws




Overall length
185 mm 7.2 inches
Round handle d=9 mm

TITANIUM

TITANIUM

Tying

VF 732




Tungsten carbide jaws

Bayonet
Overall length
200 mm 7.8 inches
Flat handle 12 mm wide


Tying

VF 733



Tungsten carbide jaws

VF 733.01





Overall length
210 mm 8.2 inches
Round handle d=9 mm

FORCEPS

Tying

VF 734




Tungsten carbide jaws


Overall length 185 mm
Round handle d=9 mm

Jewelers


VF 735




VF 735.01



VF 735.02



VF 735.03




Tungsten carbide jaws

Overall length
110 mm 4.3 inches
Flat handle 10 mm wide

Jewelers

VF 736




0.2x0.25 mm



Colibri Type
Pierse type tips
Overall length
110 mm 4.3 inches
Flat handle 10 mm wide

VF 737




0.4x0.4 mm



Tungsten carbide jaws

VF 737.01



0.4x0.35 mm
6 mm platform


Colibri Type
Overall length
120 mm 4.7 inches
Flat handle 10 mm wide

TITANIUM

TITANIUM

Jewelers

VF 738




0.3x0.25 mm



Tungsten carbide jaws

Overall length
121 mm 4.7 inches
Flat handle 10 mm wide

VF 739



0.25x0.25 mm




Tungsten carbide jaws

Overall length
135 mm 5.3 inches
Flat handle 10 mm wide


Jewelers

VF 739.01




0.25x0.25 mm
6 mm platform

VF 739.02



0.25x0.25 mm

VF 739.03




0.25x0.25 mm
6 mm platform


Tungsten carbide jaws

Different


VF 740




VF 740.01



VF 740.02




Tungsten carbide jaws




Anastomosis Type
Overall length
185 mm 7.2 inches
Round handle d=9 mm

Different


VF 741



VF 741.01




Tungsten carbide jaws




Anastomosis Type
Duble Action
Overall length
185 mm 7.2 inches
Round handle d=9 mm

Different


VF 742



VF 742.01




VF 742.02




Holding Of Fallopiev Tubes
Overall length
200 mm 7.8 inches
Round handle d=9 mm

Different


VF 743



VF 743.01



VF 743.02




Tungsten carbide jaws


Anastomosis Type
Overall length
210 mm 8.2 inches
Round handle d=9 mm

Different


VF 744




VF 744.01



VF 744.02




Tungsten carbide jaws





Anastomosis
Double Action
Overall length
210 mm 8.2 inches
Round handle d=9 mm

Different

VF 745



VF 745.01



Bayonet Spooned
Overall length
220 mm 8.6 inches
Flat handle 12 mm wide
Swilling Bayonet

TITANIUM

TITANIUM

TITANIUM



d. ext. 7.0 mm

Swelling, Bayonet
Overall length
200 mm 7.8 inches
Flat handle 12 mm wide

[illegible]


SCISSORS




Some Scissors have titanium handles

SCISSORS


VS 800




VS 800.01



VS 800.02




VS 800.03




Overall length
150 mm 5.9 inches
Round handle d=9 mm


VS 801




VS 801.01




VS 801.02





VS 801.03



VS 801.04



VS 801.05





Overall length
160 mm 6.2 inches
Flat handle 11 mm wide

STAINLESS


STAINLESS

VS 802




Double Action
Overall length
155 mm 6.1 inches
Round handle d=9 mm


VS 802.01




VS 802.02




VS 802.03





VS 802.04




VS 802.05




VS 803




VS 803.01




VS 803.02



VS 803.03




VS 803.04





Coronary scissors
Blades 14x15 mm
ball tip 1.0x0.8 mm
Overall length
175 mm 6.8 inches
Flat handle 9 mm wide

SCISSORS

VS 804





VS 804.01




Overall length
180 mm 7.0 inches
Flat handle 11 mm wide

VS 805




Stilet Type
Overall length
180 mm 7.0 inches
Round handle d=9 mm


VS 805.01




VS 805.02




VS 805.03



VS 805.04




VS 805.05





STAINLESS

TITANIUM/STAINLESS

VS 806




VS 806.01





Double Action
Left
Overall length
185 mm 7.2 inches
Round handle d=9 mm

VS 807




VS 807.01




Double Action
Right
Overall length
185 mm 7.8 inches
Round handle d=9 mm

VS 808




VS 808.01



Overall length
200 mm 7.8 inches
Flat handle 9 mm wide

VS 809




Overall length
200 mm 7.8 inches
Flat handle 11 mm wide

STAINLESS


STAINLESS

VS 810




0.3x0.1 mm

VS 810.01




0.3x0.1 mm




Double Action
Overall length
200 mm 7.8 inches
Round handle d=9 mm

VS 811



0.6x0.3 mm

VS 811.01




0.6x0.3 mm





Overall length
200 mm 7.8 inches
Round handle d=9 mm

SCISSORS

VS 812




VS 812.01




Bayonet
Overall length
200 mm 7.8 inches
Flat handle 12 mm wide

VS 813



VS 813.01




Bayonet
Overall length
200 mm 7.8 inches
Flat handle 12 mm wide


STAINLESS

STAINLESS


VS 814




VS 814.01




VS 814.02



VS 814.03




VS 814.04




Coronary
Blades 15x16 mm
tips 0.4x0.4 mm
ball tip d= 1.0x0.45 mm
Overall length
200 mm 7.8 inches
Round handle d=9 mm

VS 815



VS 815.01



Double Action
Right
Overall length
210 mm 8.2 inches
Round handle d=9 mm


STAINLESS

TITANIUM/STAINLESS


SCISSORS

TITANIUM/STAINLESS

VS 816




VS 816.01



Double Action
Left
Overall length
210 mm 8.2 inches
Round handle d=9 mm





VS 817



STAINLESS



Bayonet
Overall length
220 mm 8.6 inches
Flat handle 12 mm wide



Bayonet
Overall length
230 mm 9.0 inches
Flat handle 12 mm wide


[illegible]

VARIOUS INSTRUMENTS




TITANIUM/STAINLESS


VE 900




VE 900.01



VE 900.02




VE 900.03




Hook
Round handle d=6 mm
Overall length 200 mm


VE 901




VE 901.01




VE 901.02




VE 901.03




VE 901.04



VE 901.05




VE 901.06




Hook
Length 8.2 mm
Round handle d=6 mm
Overall length 180 mm

TITANIUM/STAINLESS


VE 902



VE 902.01



VE 902.02





Hook-retractor
Angled 90 degree., Length 25 mm
Round handle d=6 mm
Overall length 210 mm

TITANIUM/TEFLON



VARIOUS INSTRUMENTS

PR 910




Probe-retractor for gynecology
Tip 1.2 mm dia.
Round handle d=10 mm
Overall length 150 mm

PR 911



Probe-retractor for gynecology
Tip 1.2 mm dia.
Round handle d=10 mm
Overall length 185 mm

PR 912




Probe-retractor for gynecology
Tip 1.2 mm dia.
Round handle d=10 mm
Overall length 200 mm

TEFLON

TEFLON


TEFLON

PT 920




Probe for fallopiev tube
Tip 1.5 mm dia.
Overall length 300 mm


PD 930




PD 930.01




PD 930.02




PD 930.03




PD 930.04




PD 930.05




PD 930.06



PD 930.07



PD 930.08




Probes set
for
duodenum papilla


PD 931




PD 931.01



PD 931.02



PD 931.03



Bougie Vascular
Round handle d=6 mm
Overall length 180 mm


COPPER

TEFLON

STAINLESS

VARIOUS INSTRUMENTS


VL 940



Loop 4.5x2.5 mm
Round handle d=6 mm
Overall length 200 mm

TITANIUM/STAINLESS


VL 941



Loop 8.0x5.0 mm
Round handle d=6 mm
Overall length 200 mm

TITANIUM/STAINLESS


VD 950



Straight dissector
Working tip 3 mm width
Overall length 220 mm

TITANIUM/STAINLESS

VD 951





Bayonet dissector
Working tip 2 mm wide
Overall length 220 mm

CA 961



Aspiration
and
irrigation cannula
Overall length 200 mm

CA 962



Probe-cannula for introduction
contrast substance
Tip 2.0 dia.
Round handle d=8 mm
Overall length 200 mm

TITANIUM/STAINLESS

STAINLESS

STAINLESS

VARIOUS INSTRUMENTS

ZT 970

STAINLESS



Vascular clip
little
Overall length
16 mm

ZT 970.01

STAINLESS



Vascular clip
middle
Overall length
18 mm

ZT 970.02


STAINLESS



Vascular clip
big
Overall length
20 mm

VZ 980

TITANIUM



Bulldog clamp
Vascular
8 mm

VZ 980.01

TITANIUM



Bulldog clamp
Vascular
12 mm

VZ 980.02

TITANIUM



Bulldog clamp
Vascular
17 mm

RE 990



Little microsurgical retractor
Overall length 96 mm

TITANIUM

RE 991



Retractor for duodenum
Overall length 80 mm

STAINLESS

Catalogue#	Pages	Catalogue#	Pages	Catalogue#	Pages	Catalogue#	Pages
CA 961.....	54	VF 717.....	22	VH 600.04.....	3	VH 619.01.....	12
CA 962.....	54	VF 718.....	23	VH 601.....	3	VH 619.02.....	12
PD 930.....	52	VF 718.01.....	23	VH 601.01.....	3	VH 619.03.....	12
PD 930.01.....	52	VF 719.....	23	VH 602.....	4	VH 619.04.....	12
PD 930.02.....	52	VF 720.....	24	VH 602.01.....	4	VH 619.05.....	12
PD 930.03.....	52	VF 721.....	24	VH 603.....	4	VL 940.....	53
PD 930.04.....	52	VF 722.....	25	VH 603.01.....	4	VL 941.....	53
PD 930.05.....	52	VF 723.....	25	VH 604.....	5	VS 800.....	39
PD 930.06.....	52	VF 723.01.....	25	VH 604.01.....	5	VS 800.01.....	39
PD 930.07.....	52	VF 723.02.....	25	VH 604.02.....	5	VS 800.02.....	39
PD 930.08.....	52	VF 723.03.....	25	VH 604.03.....	5	VS 800.03.....	39
PD 931.....	52	VF 724.....	26	VH 605.....	5	VS 801.....	39
PD 931.01.....	52	VF 725.....	26	VH 605.01.....	5	VS 801.01.....	39
PD 931.02.....	52	VF 725.01.....	26	VH 606.....	6	VS 801.02.....	39
PD 931.03.....	52	VF 725.02.....	26	VH 606.01.....	6	VS 801.03.....	39
PR 910.....	51	VF 725.03.....	26	VH 606.02.....	6	VS 801.04.....	39
PR 911.....	51	VF 726.....	27	VH 606.03.....	6	VS 801.05.....	39
PR 912.....	51	VF 726.01.....	27	VH 606.04.....	6	VS 802.....	40
PT 920.....	52	VF 727.....	27	VH 606.05.....	6	VS 802.01.....	40
RE 990.....	55	VF 727.01.....	27	VH 607.....	6	VS 802.02.....	40
RE 991.....	55	VF 727.02.....	27	VH 607.01.....	6	VS 802.03.....	40
VD 950.....	53	VF 727.03.....	27	VH 607.02.....	6	VS 802.04.....	40
VD 951.....	54	VF 728.....	28	VH 607.03.....	6	VS 802.05.....	40
VE 900.....	50	VF 728.01.....	28	VH 607.04.....	6	VS 803.....	40
VE 900.01.....	50	VF 728.02.....	28	VH 607.05.....	6	VS 803.01.....	40
VE 900.02.....	50	VF 728.03.....	28	VH 608.....	7	VS 803.02.....	40
VE 900.03.....	50	VF 728.04.....	28	VH 608.01.....	7	VS 803.03.....	40
VE 901.....	50	VF 729.....	28	VH 609.....	7	VS 803.04.....	40
VE 901.01.....	50	VF 729.01.....	28	VH 609.01.....	7	VS 804.....	41
VE 901.02.....	50	VF 729.02.....	28	VH 609.02.....	7	VS 804.01.....	41
VE 901.03.....	50	VF 729.03.....	28	VH 609.03.....	7	VS 805.....	41
VE 901.04.....	50	VF 730.....	29	VH 609.04.....	7	VS 805.01.....	41
VE 901.05.....	50	VF 730.01.....	29	VH 609.05.....	7	VS 805.02.....	41
VE 901.06.....	50	VF 731.....	29	VH 609.06.....	7	VS 805.03.....	41
VE 902.....	50	VF 731.01.....	29	VH 609.07.....	7	VS 805.04.....	41
VE 902.01.....	50	VF 732.....	30	VH 610.....	8	VS 805.05.....	41
VE 902.02.....	50	VF 733.....	30	VH 610.01.....	8	VS 806.....	42
VF 700.....	14	VF 733.01.....	30	VH 610.02.....	8	VS 806.01.....	42
VF 700.01.....	14	VF 734.....	31	VH 610.03.....	8	VS 807.....	42
VF 701.....	14	VF 735.....	31	VH 610.04.....	8	VS 807.01.....	42
VF 701.01.....	14	VF 735.01.....	31	VH 610.05.....	8	VS 808.....	43
VF 701.02.....	14	VF 735.02.....	31	VH 611.....	8	VS 808.01.....	43
VF 701.03.....	14	VF 735.03.....	31	VH 611.01.....	8	VS 809.....	43
VF 702.....	15	VF 736.....	32	VH 612.....	9	VS 810.....	44
VF 703.....	15	VF 737.....	32	VH 612.01.....	9	VS 810.01.....	44
VF 703.01.....	15	VF 737.01.....	32	VH 613.....	9	VS 811.....	44
VF 704.....	16	VF 738.....	33	VH 613.01.....	9	VS 811.01.....	44
VF 704.01.....	16	VF 739.....	33	VH 613.02.....	9	VS 812.....	45
VF 704.02.....	16	VF 739.01.....	33	VH 613.03.....	9	VS 812.01.....	45
VF 704.03.....	16	VF 739.02.....	33	VH 614.....	10	VS 813.....	45
VF 705.....	16	VF 739.03.....	33	VH 614.01.....	10	VS 813.01.....	45
VF 706.....	17	VF 740.....	34	VH 614.02.....	10	VS 814.....	46
VF 706.01.....	17	VF 740.01.....	34	VH 615.....	10	VS 814.01.....	46
VF 706.02.....	17	VF 740.02.....	34	VH 615.01.....	10	VS 814.02.....	46
VF 707.....	17	VF 741.....	34	VH 615.02.....	10	VS 814.03.....	46
VF 708.....	18	VF 741.01.....	34	VH 615.03.....	10	VS 814.04.....	46
VF 708.01.....	18	VF 742.....	35	VH 616.....	11	VS 815.....	46
VF 709.....	18	VF 742.01.....	35	VH 616.01.....	11	VS 815.01.....	46
VF 710.....	19	VF 742.02.....	35	VH 616.02.....	11	VS 816.....	47
VF 711.....	19	VF 743.....	35	VH 616.03.....	11	VS 816.01.....	47
VF 711.01.....	19	VF 743.01.....	35	VH 616.04.....	11	VS 817.....	47
VF 712.....	20	VF 743.02.....	35	VH 616.05.....	11	VS 818.....	48
VF 713.....	20	VF 744.....	36	VH 617.....	11	VS 818.01.....	48
VF 714.....	21	VF 744.01.....	36	VH 617.01.....	11	VS 818.02.....	48
VF 715.....	21	VF 744.02.....	36	VH 617.02.....	11	VS 818.03.....	48
VF 715.01.....	21	VF 745.....	36	VH 617.03.....	11	VZ 980.....	55
VF 715.02.....	21	VF 745.01.....	36	VH 617.04.....	11	VZ 980.01.....	55
VF 715.03.....	21	VF 746.....	37	VH 617.05.....	11	VZ 980.02.....	55
VF 715.04.....	21	VH 600.....	3	VH 618.....	12	ZT 970.....	55
VF 716.....	22	VH 600.01.....	3	VH 618.01.....	12	ZT 970.01.....	55
VF 716.01.....	22	VH 600.02.....	3	VH 618.02.....	12	ZT 970.02.....	55
VF 716.02.....	22	VH 600.03.....	3	VH 619.....	12		